


Mobile X-Ray DR or Hybrid

Specification


Specifications	High Frequency X-ray Generators		
Model	POX-100BT	POXBILE-100BT	POX-100DR
Max. Output(kW)	5	5	5
kV range per 1kV step	40-110	40-110	40-110
mA Range/No. of step	20-100/6	20-100/6	20-100/6
mAs Range/No. of step	0.1-100/38	0.1-100/38	0.1-100/38
Tube Focal Spot(mm)	1.8 x 1.8	1.8 x 1.8	1.8 x 1.8
Target Angle	15°	15°	15°
SID Scale	2m	2m	2m
Battery charger	12V/24A	12V/24A	-
No. of APR	30	30	Depending on S/W

Demension


POX-100BT



POXBILE-100BT



POX-100DR



DR or Hybrid

Mobile X-Ray



Head office / R&D center / IIsan plant
 Room 405., Unitech Ville, 1141-2, Baekseok-dong, IIsandong-gu,
 Goyang, Korea 411-722
 Tel. +82-31-906-9007 (English 100~103), Fax. +82-31-908-4208
 E-mail: sales@poskom.com, www.poskom.com



Hybrid Mobile X-ray System

POX-100BT/POXBILE-100BT

POSKOM 5.0kW mobile X-ray system has powerful strength of 100mA and ideological design for mobile application. Experience high imaging resolution, user-friendly operation and smooth movement.



100mA, 5.0kW Battery Powered Mobile X-ray system for Medical

Mobile DR Solution

POX-100DR



Faster Image Solution
Higher Image Acquisition
In-House Mobility for Patient



Save, Faster therapeutic intervention

